

Tubigrip®

Multi-purpose tubular support bandage

Provides firm support for
sprains, strains and swelling



What size to choose?

Tubigrip size	Product Code		Limb circumference (cm)	Sub-bandage pressure (mmHg)
	Natural	Beige		
A	1435		10 - 14.5	20 - 30
B	1436	1449	14.5 - 29	10 - 20
C	1443	1450	22 - 35	
D	1437	1451	26 - 40	
E	1434	1448	29 - 46	
F	1438	1452	35 - 53	
G	1439	1453	40 - 62	5 - 15
J	1440		54 - 82	
K	1441		70 - 113	
L	1442		106 - 149	
M	1444		119 - 167	

Tubigrip is available under a full range of sizes to suit most applications. The sub-bandage pressures are obtained in double layer*.

Reference: Mölnlycke Health Care. Data on file.

www.molnlycke.com.au | www.molnlycke.co.nz

You will be required to input the following code **1FtYrV** or your **AHPRA registration number** to access our websites.

Mölnlycke, 12 Narabang Way, Level 4, Belrose, NSW, 2085. Phone 1800 005 231. New Zealand Orders & Enquiries 0800 005 231, www.molnlycke.co.nz. The Mölnlycke and Tubigrip trademarks, names and logos are registered globally to one or more of the Mölnlycke Health Care Group of Companies. ©2021 Mölnlycke Health Care. All rights reserved. AUWC_IM002264_V2

