

# MÖLNLYCKE - WOUND CARE PRODUCTS

PRODUCT DESCRIPTION

AREAS OF USE

SIZE

CODE

PCS/  
BOX

## BORDERED FOAM DRESSINGS



### Mepilex® Border Flex/Oval

Absorbent self-adherent foam dressing with 360° stretch Flex Technology

SafetaC  
TECHNOLOGY



- Moderate to heavily exuding wounds
- Oval shape for challenging body contours
- Diabetic foot ulcers
- Leg and foot ulcers
- Traumatic wounds, such as skin tears
- Pressure injury prevention\* and treatment

\* When a repositioning/offloading protocol is used

#### Square/Rectangle

7.5cm x 7.5cm	595211	10
10cm x 10cm	595311	10
12.5cm x 12.5cm	595011	10
15cm x 15cm	595411	10
15cm x 20cm	595611	10

#### Oval

7.8cm x 10cm	583500	5
13cm x 16cm	583300	5
15cm x 19cm	583400	5



### Mepilex® Border Heel

Absorbent self-adherent foam dressing for the heel area

#### Added benefits:

- Thicker borders to minimise rolling and to allow efficient handling
- Optimised seal for enhanced protection of the high-risk heel area
- Additional handling tabs for easier skin checks
- Enhanced coverage of high risk areas

SafetaC  
TECHNOLOGY



- Moderate to heavily exuding wounds for the heel area
- Heel pressure injury prevention\* and treatment (all stages)
- Diabetic foot ulcers
- Leg ulcers
- Traumatic wounds

\* When a repositioning/offloading protocol is used

22cm x 23cm	282750	6
-------------	--------	---



### Mepilex® Border Sacrum

Absorbent self-adherent foam dressing for the sacral area

#### Added benefits:

- Thicker borders to minimise rolling and to allow efficient handling
- Optimised gluteal seal for enhanced protection of the high-risk sacral area
- Enhanced coverage of high risk areas
- Additional handling tabs for easier skin checks


SafetaC  
TECHNOLOGY



- Moderate to heavily exuding wounds for the sacral area
- Sacral pressure injury prevention\* and treatment (all stages)

\* When a repositioning/offloading protocol is used

16cm x 20cm	282050	5
22cm x 25cm	282450	5



### Mepilex® Border Flex Lite

Thin absorbent self-adherent foam dressing with 360° stretch Flex Technology

SafetaC  
TECHNOLOGY



- Non- to moderately exuding wounds
- Leg and foot ulcers
- Pressure injuries
- Surgical and traumatic wounds, e.g. abrasions, blisters and skin tears
- Protection of compromised and/or fragile skin

4cm x 5cm	581011	10
5cm x 12.5cm	581100	5
7.5cm x 7.5cm	581200	5
10cm x 10cm	581300	5
15cm x 15cm	581500	5

## REVOLUTIONISE CARE FOR PEOPLE AND PLANET.



Visit [my.molnlycke.com.au](https://my.molnlycke.com.au) or scan the QR code for product information, user instructions and more

  
**Mölnlycke®**

# MÖLNLYCKE - WOUND CARE PRODUCTS

PRODUCT DESCRIPTION


AREAS OF USE

SIZE

CODE

PCS/  
BOX

## NON-BORDERED FOAM DRESSINGS



### Mepilex® Up

Innovative absorbent foam dressing, with a patented foam structure specifically designed to minimise the risk of leakage

**SafetaC**  
TECHNOLOGY



- Low to heavily exuding wounds in all healing stages
- Hard-to-heal (Chronic) wounds
- Leg ulcers and foot ulcers
- Pressure injuries
- Traumatic wounds, e.g. skin tears
- Can be used under compression

5cm x 5cm	212015	5
10cm x 10cm	212100	5
15cm x 15cm	212300	5
10cm x 20cm	212200	5
20cm x 20cm	212400	5
20cm x 50cm	212500	2



### Mepilex®

Highly conformable and absorbent foam dressing

**SafetaC**  
TECHNOLOGY



- Low to moderately exuding wounds
- Leg ulcers and diabetic foot ulcers
- Pressure injuries
- Traumatic wounds
- Partial-thickness burns
- Can be used under compression

5cm x 5cm	294015	5
10cm x 10cm	294100	5
10cm x 20cm	294200	5
15cm x 15cm	294300	5
20cm x 20cm	294400	5
20cm x 50cm	294500	2



### Mepilex® Heel

Absorbent foam dressing for the heel area

**SafetaC**  
TECHNOLOGY



- A wide range of exuding wounds on the heel or malleolus
- Foot ulcers / diabetic foot ulcers

13cm x 20cm	288100	5
-------------	--------	---



### Mepilex® Lite

Thin absorbent foam dressing

**SafetaC**  
TECHNOLOGY




- Non- to low exuding wounds
- Pressure injuries
- Diabetic foot ulcers
- Leg ulcers
- Radiation skin reactions
- Epidermolysis Bullosa

6cm x 8.5cm	284000	5
10cm x 10cm	284100	5
15cm x 15cm	284300	5
20cm x 50cm	284500	4
12.5cm x 12.5cm*	284123	5

\* 12.5cm x 12.5cm cutting guide

## WOUND CONTACT LAYERS



### Mepitel® One

Wound contact layer dressing with **SafetaC**® on one side

**SafetaC**  
TECHNOLOGY



- Skin tears
- Lacerations
- Surgical incisions
- Fixation of partial- and full-thickness grafts
- Partial-thickness burns
- Diabetic foot ulcers
- Venous and arterial leg ulcers

5cm x 7.5cm	289100	10
6cm x 7cm*	289170	5
7.5cm x 10cm	289300	10
9cm x 10cm*	289270	5
13cm x 15cm*	289470	5
10cm x 18cm	289500	10
17cm x 25cm	289700	5
24cm x 27.5cm	289670	5
27.5cm x 50cm	289750	5

\* Sizes with pre-printed symbols



### Mepitel®

Wound contact layer dressing with **SafetaC**® on both sides

**SafetaC**  
TECHNOLOGY



- Skin tears
- Surgical incisions
- Partial-thickness burns
- Fixation of partial- and full-thickness grafts
- Diabetic foot ulcers, venous and arterial leg ulcers

5cm x 7.5cm	290510	10
7.5cm x 10cm	290710	10
10cm x 18cm	291010	10
20cm x 30cm	292005	5



### Mepilex® Transfer

Exudate transfer foam dressing

**SafetaC**  
TECHNOLOGY



- Highly exuding wounds
- Difficult-to-dress areas
- Radiation skin reactions
- Epidermolysis Bullosa

7.5cm x 8.5cm	294600	5
10cm x 12cm	294700	5
15cm x 20cm	294800	5
20cm x 50cm	294502	4

# MÖLNLYCKE - WOUND CARE PRODUCTS

PRODUCT DESCRIPTION

AREAS OF USE

SIZE

CODE

PCS/  
BOX

## POST-OPERATIVE DRESSINGS



### Mepilex® Border Post-Op

Self-adherent post-operative dressing with 360° stretch Flex Technology

**SafetaC**  
TECHNOLOGY



- Low to heavily exuding wounds
- Optimised for post-operative use and blood absorption
- Surgical wounds, cuts and abrasions
- Ideal for use over joints, such as hips and knees

6cm x 8cm	496100	10
9cm x 10cm	496200	10
10cm x 15cm	496300	10
10cm x 20cm	496400	10
10cm x 25cm	496450	10
10cm x 30cm	496600	10
10cm x 35cm	496650	5

## FILM DRESSINGS



### Mepitel® Film

Breathable, transparent film dressing with SafetaC® Technology

**SafetaC**  
TECHNOLOGY

- Protects fragile and sensitive skin
- For fixation of primary dressings
- For use under aggressive tapes
- For use in radiotherapy
- Can be used as a protective dressing over IV sites

6cm x 7cm	296170	10
10cm x 12cm	296270	10
10cm x 25cm	296470	10
15cm x 20cm	296670	10



### Mepore® Film

Breathable, transparent film dressing

- Designed as a secondary fixation dressing and as a primary dressing for a wide range of clean wounds

6cm x 7cm	270600	100
10cm x 12cm	271500	70
10cm x 25cm	272500	20
15cm x 20cm	273000	10

## ISLAND DRESSINGS



### Mepore®

Absorbent island dressing



- For a wide variety of superficial wounds which have low to moderate levels of exudate
- Post-operative wounds
- Cuts and abrasions

5cm x 7cm	672000	60
6cm x 7cm	670800	60
8cm x 10cm	672100	50
9cm x 10cm	670900	50
8cm x 15cm	672200	50
9cm x 15cm	671000	50
9cm x 20cm	671100	30
9cm x 25cm	671200	30
9cm x 30cm	671300	30
9cm x 35cm	671400	30



### Mepore® Pro

Shower-proof, absorbent island dressing



- For a wide variety of superficial wounds which have low to moderate levels of exudate
- Where a shower-proof dressing is required
- Post-operative wounds
- Cuts and abrasions

6cm x 7cm	670820	60
9cm x 10cm	670920	40
9cm x 15cm	671020	40
9cm x 20cm	671120	30
9cm x 25cm	671220	30
9cm x 30cm	671320	30



### Mepore® Film&Pad

Transparent, shower-proof absorbent island dressing



- For a wide variety of superficial wounds which have low to moderate exudate levels
- Where a shower-proof dressing is required
- Post-operative wounds
- Cuts and abrasions

4cm x 5cm	275100	85
5cm x 7cm	275200	85
<b>Oval</b>		
5cm x 7cm	275300	85
<b>Square</b>		
9cm x 10cm	275400	30
9cm x 15cm	275500	30
9cm x 20cm	275600	30
9cm x 25cm	275700	30
9cm x 35cm	275900	25

# MÖLNLYCKE - WOUND CARE PRODUCTS

PRODUCT DESCRIPTION

AREAS OF USE

SIZE

CODE

PCS/  
BOX

## ANTIMICROBIAL DRESSINGS



### Mepilex® Border Ag

Antimicrobial absorbent self-adherent foam dressing

*SafetaC*  
TECHNOLOGY



- Moderate to heavily exuding wounds
- Infected wounds
- Leg and foot ulcers
- Pressure injuries
- Malignant wounds
- Surgical incisions

7.5cm x 7.5cm

395200

5

10cm x 10cm

395300

5

15cm x 15cm

395400

5

15cm x 20cm

395600

5

**Surgical Sizing**

10cm x 20cm

395800

5

10cm x 25cm

395700

5

10cm x 30cm

395900

5



### Mepilex® Border Sacrum Ag

Antimicrobial absorbent self-adherent foam dressing for the sacral area

*SafetaC*  
TECHNOLOGY



- Moderate to heavily exuding wounds for the sacral area where a self-adherent antimicrobial dressing is required

18cm x 18cm

382000

5

23cm x 23cm

382400

5



### Mepilex® Ag

Antimicrobial absorbent foam dressing for the sacral area

*SafetaC*  
TECHNOLOGY



- Low to moderately exuding wounds
- Partial-thickness burns
- Infected wounds
- Leg ulcers and foot ulcers
- Pressure injuries
- Can be used under compression

10cm x 10cm

287110

5

10cm x 20cm

287210

5

15cm x 15cm

287310

5

20cm x 20cm

287410

5

20cm x 50cm

287510

2



### Mepilex® Heel Ag

Antimicrobial absorbent foam dressing for the heel area

*SafetaC*  
TECHNOLOGY




- Low to moderately exuding wounds
- Infected wounds on the foot or heel
- Diabetic foot ulcers
- Pressure injuries
- Can be used under compression

13cm x 20cm

388100

5



### Mepilex® Transfer Ag

Antimicrobial compressed foam, exudate transfer dressing

*SafetaC*  
TECHNOLOGY



- Heavily exuding wounds
- Infected or malignant wounds
- Donor sites
- Partial-thickness burns
- Leg ulcers and diabetic foot ulcers
- Difficult-to-dress areas

7.5cm x 8.5cm

394000

10

10cm x 12.5cm

394100

5

12.5cm x 12.5cm

394700

5

15cm x 20cm

394800

10

20cm x 50cm

394500

2



### Melgisorb® Ag

Highly absorbent antimicrobial alginate dressing



- Moderate to heavily exuding wounds
- Post-operative wounds
- Leg ulcers and diabetic foot ulcers
- Trauma wounds (dermal lesions, trauma injuries or incisions)
- Graft and donor sites
- Pressure injuries

5cm x 5cm

256055

10

10cm x 10cm

256105

10

15cm x 15cm

256155

10

20cm x 30cm

256205

5

3cm x 44cm

256605

10

(Cavity)



### Exufiber® Ag+

Antimicrobial gelling fibre dressing with Hydrolock® Technology



In the following moderate to high exuding wounds:

- Venous leg ulcers
- Diabetic foot ulcers
- Pressure injuries

5cm x 5cm

603401

10

10cm x 10cm

603402

10

15cm x 15cm

603403

10

20cm x 30cm

603407

5

4.5cm x 20cm

603405

10

2cm x 45cm

603400

5

(Cavity)

# MÖLNLYCKE - WOUND CARE PRODUCTS

PRODUCT DESCRIPTION

AREAS OF USE

SIZE

CODE

PCS/  
BOX

## SUPERABSORBENTS



### Mextra® Superabsorbent

Superabsorbent dressing with a fluid-repellent backing



- Moderate to very heavily exuding wounds
- Can be used as a primary or secondary dressing where high absorbency is required, e.g. together with Mepitel®, Mepitel One or Mepilex® Transfer
- Use with Tubifast® for light retention

12.5cm x 12.5cm	610000	10
12.5cm x 17.5cm	610100	10
12.5cm x 22.5cm	610200	10
17.5cm x 22.5cm	610300	10
22.5cm x 27.5cm	610400	10
22.5cm x 32.5cm	610500	10
22.5cm x 42.5cm	610600	10



### Mesorb®

Highly absorbent dressing with a fluid-repellent backing



- Moderate to very heavily exuding wounds
- Can be used as a primary or secondary dressing where high absorbency is required, e.g. together with Mepitel, Mepitel One or Mepilex Transfer
- Use with Tubifast for light retention

10cm x 13cm	677001	50
13cm x 15cm	677201	50
10cm x 23cm	677401	50
15cm x 23cm	677601	50
23cm x 25cm	677701	30
23cm x 30cm	677801	30

## FIBRE/ALGINATE DRESSINGS



### Exufiber®

Gelling fibre dressing with Hydrolock® Technology



- Leg ulcers and diabetic foot ulcers
- Pressure injuries
- Partial-thickness burns
- Surgical wounds
- Donor sites
- Malignant wounds

5cm x 5cm	709900	10
10cm x 10cm	709901	10
15cm x 15cm	709903	10
20cm x 30cm	709904	5
4.5cm x 20cm	709906	10
1cm x 45cm (Cavity)	709908	5
2cm x 45cm (Cavity)	709909	5



### Melgisorb® Plus

Highly absorbent alginate dressing



- Moderate to heavily exuding wounds
- Pressure injuries
- Venous and arterial leg ulcers
- Diabetic foot ulcers
- Post-operative wounds
- Dermal lesions and other external wounds inflicted by trauma

5cm x 5cm	252000	10
10cm x 10cm	252200	10
10cm x 20cm	252500	10
3cm x 45cm (Cavity)	253500	5

## WOUND IRRIGATION



### Granudacyn®

Wound irrigation solution and gel for cleansing and moisturising acute, hard-to-heal (chronic) and contaminated wounds

- All hard-to-heal (chronic) wounds of any depths, such as diabetic foot ulcers, pressure injuries, venous leg ulcers etc.
- All acute wounds, such as cuts, bite wounds, lacerations and abrasions
- Surgical wounds (intra-operative and post-operative)
- Wounds with exposed cartilage, tendons, ligaments and/or bones
- Superficial and partial-thickness burns
- Radiation ulcers
- Fistulas and abscesses

#### Solution

50ml spray	360150	1
250ml spray	360100	1
500ml	360101	1
1000ml	360102	1

#### NPWT Solution

500ml	360103	1
1000ml	360104	1

#### Gel

50g	360107	1
100g spray	360108	1

# MÖLNLYCKE - WOUND CARE PRODUCTS

PRODUCT DESCRIPTION

AREAS OF USE

SIZE

CODE

PCS/  
BOX

## BANDAGES



### Tubifast®

Lightweight tubular retention bandage

WITH  
2-WAY STRETCH®  
TECHNOLOGY

- Lightweight tubular dressing retention bandage that is ideal for any part of the body
- Can be used for patch wrapping and as an undercast stockinette
- Can be used under compression
- For use as skin covering
- Epidermolysis Bullosa

10m roll

RED	3.5cm x 10m	2434	1
GREEN	5cm x 10m	2436	1
BLUE	7.5cm x 10m	2438	1
YELLOW	10.75cm x 10m	2440	1
PURPLE	20cm x 10m	2444	1



### Tubifast® Garments

Ready-to-wear garments

WITH  
2-WAY STRETCH®  
TECHNOLOGY

- Ready-to-wear garments for dressing retention and wet or dry wrapping, are available in full-sleeved vests, tights, leggings, socks and gloves
- For children of six months to 14 years of age
- With no tapes, pins or ties, Tubifast Garments are quick and easy to apply, and their 2-Way Stretch construction ensures complete freedom of movement, so aiding compliance.
- Comfortable to wear under nightwear and ordinary clothes

Vests

6-24 months	992007	6
2-5 yr	992008	6
5-8 yr	992009	6
8-11 yr	992010	6
11-14 yr	992011	6

Tights

6-24 months	992012	6
-------------	--------	---

Leggings

2-5 yr	992013	6
5-8 yr	992014	6
8-11 yr	992015	6
11-14 yr	992016	6

Socks

1 size: 2-14yr	992017	6
----------------	--------	---

Gloves

Small child	5923	6
-------------	------	---



### Tubigrip®

Multi-purpose, elasticated tubular support bandage that provides firm support for sprains, strains and swelling

Size / Limb

Circumference (cm)

	Natural (N)	Beige (B)	N	B
A 10 - 14.5	1435	-	1	-
B 14.5 - 29	1436	1449	1	1
C 22 - 35	1443	1450	1	1
D 26 - 40	1437	1451	1	1
E 29 - 46	1434	1448	1	1
F 35 - 53	1438	1452	1	1
G 40 - 62	1439	1453	1	1
J 54 - 82	1440	-	1	-
K 70 - 113	1441	-	1	-
L 106 - 149	1442	-	1	-
M 119 - 167	1444	-	1	-



### Setopress®

Lightweight high compression bandage

- Used in the treatment of venous leg ulcers and associated conditions requiring graduated compression of the leg

10cm x 3.5m	3505	12
-------------	------	----



# MÖLNLYCKE - WOUND CARE PRODUCTS

PRODUCT DESCRIPTION

AREAS OF USE

SIZE

CODE

PCS/  
BOX

## FIXATION



### Mepitac®

Polyurethane tape with Safetac® Technology

Safetac®  
TECHNOLOGY



- For the fixation of medical devices, such as drains, tubes and probes
- Skin protection under devices

2cm x 3m  
4cm x 1.5m

298300  
298400

1  
1



### Mefix®

Self-adhesive fabric fixation tape

- For fixation of a wide variety of dressings, tubes and cannulas

10m roll

x 2.5cm  
x 5cm  
x 10cm  
x 15cm  
x 20cm  
x 30cm

310299  
310599  
311099  
311599  
312000  
313000

1  
1  
1  
1  
1  
1

## SCAR THERAPY & DERMATOLOGY



### Mepiform®

Scar dressing with Safetac® Technology

Safetac®  
TECHNOLOGY

- Medical grade silicone dressing for the prevention and treatment of all scars, in particular keloid scars and/or hypertrophic scars

4cm x 30cm  
5cm x 7.5cm  
10cm x 18cm

293100  
293200  
293400

5  
5  
5



### Epaderm® Cream

An emollient for dry skin conditions

### Epaderm® Junior Cream\*

- A 2-in-1 emollient that can be used as a skin emollient, as well as a skin cleanser
- Used for eczema, psoriasis and all other dry skin conditions
- Suitable for all ages with the same formulation for both adults and children

25g  
50g  
500g  
50g\*  
50g\*

99400822  
99400851  
99400850  
99400860  
99400861

12  
12  
6  
12  
6



### Epaderm® Ointment

An emollient for dry skin conditions

### Epaderm® Junior Ointment\*

- A 3-in-1 emollient that can be used as a skin emollient, a skin cleanser and a bath additive
- Used for eczema, psoriasis and all other dry skin conditions
- Suitable for all ages with the same formulation for both adults and children

25g  
125g  
500g  
1kg  
125g\*

400861  
99400846  
99400847  
99400818  
99400871

12  
6  
6  
6  
6

## ACTIVE DRESSINGS



### Mesalt®

Sodium chloride impregnated dressing



- For the management of moderate to heavily exuding non-infected/infected wounds, such as pressure injuries, foot ulcers and deep wound management, e.g. cavity wounds and fistulas

5cm x 5cm  
7.5cm x 7.5cm  
10cm x 10cm  
2cm x 1m  
(Cavity)

285580  
285780  
286080  
285280

30  
30  
30  
10

# MÖLNLYCKE - WOUND CARE PRODUCTS

PRODUCT DESCRIPTION

AREAS OF USE

SIZE

CODE

PCS/  
BOX

## TURNING & POSITIONING SYSTEMS



### Z-Flex™ Heel Boot

Assists with offloading of the heel

- Fully offloads the heel
- Positions feet in neutral alignment
- Accommodates use of compression devices
- Built-in gate allows quick access for heel inspection
- One size fits most

2 boots per box  
8 boots per box

1400122  
1400123

2  
8



### Z-Flo™ Fluidized Positioner

Assists with patient positioning and offloading of at risk anatomical sites

- Redistributes pressure or offloads
  - Keeps patients in position
  - Mouldable to protect vulnerable areas
- (More sizes available for Paediatric and Neonatal patients)

#### With Tab

40cm x 76cm

1401003

1

#### Without Tab

64cm x 91cm

1401005

1

30cm x 51cm

1401007

8

(Occipital)



### Tortoise® Turning & Positioning System

Assists with turning and repositioning of patients

- Facilitates patient turning, positioning, boosting and lateral transfers
- Helps maintain the intended degree of turn and protects bony prominences
- Reduces the risk of caregiver injury

#### Standard

96.5cm x 147cm

1400800

1

#### Bariatric

124.5cm x 147cm

1400801

1

## NEGATIVE PRESSURE WOUND THERAPY



### Avance® Solo

A portable single-patient use negative pressure (NPWT) system for closed surgical incisions and open wounds


- -125mmHg negative pressure for up to 14 days
- Visual and audible alarms for blockage, leakage and low battery
- Excess fluid transferred into a 50ml replaceable canister
- CFM® Technology ensures continuous negative pressure
- Absorbent multilayer dressing with Safetac® Technology

**Note:** all sets/kits containing a pump also include 4 x AA batteries and a belt clip.

\* This indicates complete sets that contain a pump, a canister and a dressing.

10cm x 20cm*	8881020	1
10cm x 30cm*	8881030	1
10cm x 35cm*	8881035	1
Pump & canister kit	8880050	2
10cm x 20cm	881020	2
10cm x 30cm	881030	2
10cm x 35cm	881035	2
15cm x 15cm	881515	2
15cm x 20cm	881520	2
15cm x 30cm	881530	2
20cm x 20cm	882020	2
25cm x 25cm	882525	2
50ml canister	880050	4
Avance Solo Foam	882000	4
1.5cm x 10cm x 12cm		
Carrying Case	890003	2
Carry Bag Small	890001	1

Safetac®  
TECHNOLOGY



### Avance® Solo Adapt

A portable single-patient use negative pressure (NPWT) system for pressure injuries

- -125mmHg pressure for up to 14 days
- Visual and audible alarms for blockage, leakage and low battery
- Excess fluid transferred into a 50ml replaceable canister
- Controlled Fluid Management® Technology ensures continuous negative pressure therapy
- Single use polyurethane foam wound filler which can be cut to the shape of the wound
- Avance Solo Adpat Film with Safetac® Technology

Avance Solo Adapt Set	9994152	1
Pump Cannister, Film, Foam and 4 x Batteries		
Avance Solo Adapt Dressing Kit	994152	5
Avance Solo Adapt Film	994000	5
20 x 27cm		
Carrying Case	890003	2
Carry Bag Small	890001	1

Safetac®  
TECHNOLOGY



NB: Please refer to the Instructions of Use prior to usage.

Find out more at [www.molnlycke.com.au](http://www.molnlycke.com.au)

You will be required to input the following code **1FtYrV** or your **AHPRA registration number** to access our websites.

Mölnlycke Health Care, 12 Narabang Way, Level 4, Belrose, NSW. 2085, Australia. Phone 1800 005 231. The Mölnlycke, 2-Way Stretch, Avance, Epaderm, Exufiber, Granudacyn, Hydrolock, Mefix, Melgisorb, Mepiform, Mepilex, Mepitac, Mepitel, Mepore, Mesalt, Mesorb, Mextra, Tubifast, Tubigrip, Safetac, Setopress, Tortoise names and logos are registered globally to one or more of the Mölnlycke Health Care Group of Companies. Z-Flex and Z-Flo are trademarks used by Mölnlycke Health Care Group of Companies.

© 2025 Mölnlycke Health Care. All rights reserved. AUWC\_001704

